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



laser based sensor systems

Fachtagung Bauwerksdiagnose 2018

Analytics of Minerals	Offline analytics	Inline XRF (Gamma in, Gamma out)	Inline Neutrons (Neutron in Gamma out)	Inline LIBS (Heat in Photons out)
Speed of result	Very slow	fast	slow	fast
elements	all	Light difficulty	Heavy difficulty	all
precision	Very good	Semi quantitative	qualitative	Typ. 5% rel.
Lateral resolution	meters	>5mm	meters	100µm
Measurements/min	< 1	60	1-2	6000
Safety regulations	no	high	high	Low if in housing
others	Sampling error, sampling costs	Distance sensitive	much material needed	Just point by point

# LIBS Basics


**SECOPTA**  
laser based sensor systems

Fachtagung Bauwerksdiagnose 2018

3

## Each Element has a unique fingerprint

**SECOPTA**  
laser based sensor systems

Fachtagung Bauwerksdiagnose 2018

4

# LIBS Basics – Advantages of the method

## Universal and flexible

- each element is detected, even lightweight
- identification of complex objects by “fingerprint”
- adjustable by software

Quantitative

Qualitative

Process qualified

## Fast and straight forward

- **Online** – full automated data analysis in ~1 ms  
FiberLIBS with up to 100 measurements per second  
applicable for monitoring und process control
- **Inline** – compact sensor head for integration directly into machinery
- **Insitu** – no sampling and no sample preparation

## USPs SECOPTA LIBS Technology


fast.

- Results within milliseconds
- > 350 measurements / s
- Up to 20.000 measurement points / s

precise.

- Quantitative Measurements <0,1 wt%
- Typical error of approx. 5% rel. possible @ 3m/s
- Simultaneous measurement of all elements (even lightweight)

robust.

- SECOPTA sensors in industry quality
- IP67, sealed in clean room environment, air stream
- Connection to PLC via common industrial interfaces



# LIBS solutions from SECOPTA: Even more USPs

- Non-contact-measuring
- Product feature: autofocus
- Works under extreme conditions (e.g. directly on liquid metal)
- Directly on conveyor belt / in motion
- No sample preparation
- Element mapping: additional source of information
- Continuous , **100% quality control** (diversification!)
- Cleaning and measuring may be one production step (PA feature)
- Ready for industry 4.0
- Team of experts
- SECOPTA team and Japanese Partners provides exceptional service
- Consistent product quality Made in Germany

## Quantitative calibration

- Quantitative calibration
  - ✓ Multivariate: Use the hole spectrum, PLS (Partial least square regression)
  - ✓ Pattern recognition in complex data (fingerprints)
- Advantages over conventional methods
  - ✓ Using all lines – fingerprint
  - ✓ Solving problems with matrix effects

# LIBS solutions from Secopta

**SECOPTA**  
laser based sensor systems

Fachtagung Bauwerksdiagnose 2018

9

## Comparision Spectrometers for LIBS


Spectrometers for LIBS	Echelle + ICCD	Paschen-Runge	Compact (Tscherney Turner)
Speed	20 Hz	1 kHz	800 Hz
resolution	< 20 pm	50 pm	0.15 nm
Stability under industrial conditions	poor	ok	good
Price	high	mid	low
others	High sensitive	Combination of photomultipliers and CCD lines	easy to stabilize temperature

**SECOPTA**  
laser based sensor systems

Fachtagung Bauwerksdiagnose 2018

Laser for LIBS	High Power (NdYAG)	Low Power / short pulse	FiberLaser (1ns)	FiberLaser (100ns) + Sparc
Repetition rate	80Hz	100Hz	10kHz	<b>20kHz</b>
Pulse energy	25-1000mJ	3mJ	1mJ	1mJ
Beam quality	good	Very good	good	good
Size/weight	big / >10kg	3kg	4kg	2 kg
Long Term stability	< 1year	1-3 years	?	<b>10 years</b>
Price	high	middle	high	<b>low</b>
others	Degeneration while working	Not stable within first shots	Not a standard product yet	<b>Well establish in laser marking</b>

## FiberLIBS lab – MK II

- Offline material analysis
- Detection of contaminations
- Determination of layer thickness 1D – 2D – 3D Scans
- Element-Mapping
- Depth- profiles

### Features

Laser Distance Measurement  
Full automated 3D-Scanning


Software Package: Chemometrics  
2D-Elemental Mapping

No problems with Laser Safety

## Example: 2D- Element mapping

**Analysis of Fe and Si concentration  
depending on mineralic composition of the sample**

## Summary

-  LIBS is an innovative universal elemental analysis technique
-  fast. LIBS is extreme fast, results within ms, >350 measurements/s
-  precise. LIBS is much more precise than other process measurement techniques like XRF or neutrons
-  robust. LIBS systems from Secopta are extremely robust: Low maintenance
-  LIBS can be used under harsh industrial environmental conditions, e.g. for primary raw materials - sorting or quality control
-  LIBS can be used quantitative or qualitative (classification), Secopta provides standard applications as well as individual chemometric models

## Customers Benefit

High add Value for the Customer  
Low Cost of Ownership

Reliable long term Partnership is our Key Value

Our local Partners are essential for our sales approach

Professional Service together with our  
Japanese - German Service Team

Easy and mediate Remote maintenance possible upon request

High Quality - Made in Germany

**SECOPTA**

laser based sensor systems

Fachtagung Bauwerksdiagnose 2018

15

Thank You! Contact us...

**SECOPTA**  
laser based sensor systems

Fachtagung Bauwerksdiagnose 2018

16